

ESAW Medium

(Enriched Seawater, Artificial Water)

This is a medium designed to mimic natural seawater that coastal and open ocean phytoplankton live in (e.g. diatoms *Thalassiosira pseudonana* aka "THAPS").

This recipe is for a final volume of 1L medium (you may multiple the recipe for any quantity you need)

- Determine the final volume of media you want to create and calculate the proper grams and volumes to add.
- Collect all needed materials:
 - 2 mixing jars (or flasks) large enough to split the final volume of media in (eg. If making 1L of media, use 1L flask and a 500mL flask)
 - graduated cylinder filled with DI H2O to proper final volume
 - weigh boats and spatulas
 - gloves (optional)
 - Erlenmeyer flask (2L) that is sterile
 - Filter tops to filter media through (or if you do not have filters, see other techniques in the media prep folder)
 - labeling tape
 - pen
- Dissolve Salt I (anhydrous salts) first in ~600mL of DI water
 - Weigh out each salt and place in flask (don't forget to tare the scale between each weigh)
 - **Before adding the stock salts, dissolve dry salts in the water.

Once dissolved, then add the stock salt solutions.

Component	Stock Solution	Quantity in 1L final medium	Molar Concentration in Final
Salt I			
Anhydrous salts			
NaCl		21.194 g	3.63×10^{-1}
Na ₂ SO ₄		3.550 g	2.50×10^{-2}
KCl		0.599 g	8.03×10^{-3}
NaHCO ₃		0.174 g	2.07×10^{-3}
KBr**	stock	**1mL	7.25×10^{-4} M
H ₃ BO ₃	stock	**1mL	3.72×10^{-4} M
NaF	stock	**1mL	6.67×10^{-3} M

86.3g/L
23g/L
2.8g/L

- Dissolve Salt II (hydrated salts) in ~400mL of water

Salt II	Stock Solution	Quantity	Molar Concentration in Final
Hydrated salts			
MgCl ₂ 6H ₂ O		9.592 g	4.71×10^{-2} M
CaCl ₂ 2H ₂ O		1.344 g	9.14×10^{-3} M
SrCl ₂ 6H ₂ O	stock	**1mL	8.18×10^{-5} M

21.89g/L

- Combine dissolved salt I and II, and then add Trace Metals, Vitamins and Major Nutrients (NO₃, PO₄, SiO₃).

**Pay attention to whether you are creating a media that is limited in one nutrient!!

Trace Metals	1mL
Vitamins	1mL

For full "f/2" Media				For Nitrate limiting Media				For Silica limiting Media				For Phosphate limiting Media			
Major Nutrients				Major Nutrients				Major Nutrients				Major Nutrients			
NaNO ₃	75 g L ⁻¹ dH ₂ O	1 mL	88uM	NaNO ₃	75 g L ⁻¹ dH ₂ O	0.5mL	44uM	NaNO ₃	75 g L ⁻¹ dH ₂ O	1 mL	88uM	NaNO ₃	75 g L ⁻¹ dH ₂ O	1 mL	88uM
NaH ₂ PO ₄ ·H ₂ O	5 g L ⁻¹ dH ₂ O	1 mL	36.2uM	NaH ₂ PO ₄ ·H ₂ O	5 g L ⁻¹ dH ₂ O	1 mL	36.2uM	NaH ₂ PO ₄ ·H ₂ O	5 g L ⁻¹ dH ₂ O	1 mL	36.2uM	NaH ₂ PO ₄ ·H ₂ O	5 g L ⁻¹ dH ₂ O	0.071mL	2.5uM
Na ₂ SiO ₃ ·9H ₂ O	30 g L ⁻¹ dH ₂ O	1mL	106uM	Na ₂ SiO ₃ ·9H ₂ O	30 g L ⁻¹ dH ₂ O	1mL	106uM	Na ₂ SiO ₃ ·9H ₂ O	30 g L ⁻¹ dH ₂ O	0.333mL	37.5uM	Na ₂ SiO ₃ ·9H ₂ O	30 g L ⁻¹ dH ₂ O	1mL	106uM

- Filter sterilization is recommended. (See Media Prep folder for more information if needed.)
- Label the media with [ESAW] [condition (f/2, N-, Si- or P-)] [date]
- Let media sit overnight before use. The final pH is 8.2.